



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

AllView SDS-PAGE Buffer ABN-TM0011

Article Name	AllView SDS-PAGE Buffer
Biozol Catalog Number	ABN-TM0011
Supplier Catalog Number	TM0011
Alternative Catalog Number	ABN-TM0011-500
Manufacturer	Abnova
Category	Sonstiges
Application	EP
Product Description	AllView SDS-PAGE Buffer has a remarkable feature that enable us to separate proteins with wide range of molecular weights in the basic Laemmli gel (Tris-HCl) similar to using gradient gel. After electrophoresis, the polyacrylamide gel can be used dir
Form	Liquid
Application Dilute	Electrophoresis (dilute with ultrapure water to 1x solution before use).