

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### 5p Subtelomere(DEAC) FISH Probe ABN-FE0044

Article Name	5p Subtelomere(DEAC) FISH Probe
Biozol Catalog Number	ABN-FE0044
Supplier Catalog Number	FE0044
Alternative Catalog Number	ABN-FE0044-20
Manufacturer	Abnova
Category	Sonstiges
Application	FISH
Species Reactivity	Human
Product Description	Labeled FISH probes for identification of subtelomere aberrations using Fluorescent In Situ Hybridization Technique. ( <a href="#">Technology</a> )....
Application Dilute	The probe is provided in 5x concentrated format, to allow mixing of up to 5 Subtelomere FISH Probes in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer ( <a href="http://www.abnova.com/products">Cat &amp; href=http://www.abnova.com/products</a>