

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

19p Subtelomere(DEAC) FISH Probe ABN-FE0169

| | |
|----------------------------|--|
| Article Name | 19p Subtelomere(DEAC) FISH Probe |
| Biozol Catalog Number | ABN-FE0169 |
| Supplier Catalog Number | FE0169 |
| Alternative Catalog Number | ABN-FE0169-20 |
| Manufacturer | Abnova |
| Category | Sonstiges |
| Application | FISH |
| Species Reactivity | Human |
| Product Description | Labeled FISH probes for identification of subtelomere aberrations using Fluorescent In Situ Hybridization Technique. (Technology).... |
| Application Dilute | The probe is provided in 5x concentrated format, to allow mixing of up to 5 Subtelomere FISH Probes in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer (Cat & a href=http://www.abnova.com/products |