



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GloMelt(TM) Thermal Shift Protein Stability Kit with ROX BOT-33022-T

Article Name	GloMelt(TM) Thermal Shift Protein Stability Kit with ROX
Biozol Catalog Number	BOT-33022-T
Supplier Catalog Number	33022-T
Alternative Catalog Number	BOT-33022-T
Manufacturer	Biotium
Category	Kits/Assays
Product Description	GloMelt(TM) dye undergoes fluorescence enhancement upon binding to hydrophobic regions of denatured proteins. GloMelt(TM) can therefore be used to detect protein unfolding or protein thermal stability by performing a thermal shift assay, also called
Source	Animal