



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **GelRed Prestain Plus 6X DNA Loading Dye BOT-41011**

| Article Name               | GelRed Prestain Plus 6X DNA Loading Dye  |
|----------------------------|--|
| Biozol Catalog Number      | BOT-41011  |
| Supplier Catalog Number    | 41011  |
| Alternative Catalog Number | BOT-41011  |
| Manufacturer               | Biotium  |
| Category                   | Molekularbiologie  |
| Product Description        | This 6X DNA loading buffer includes Biotiums non-mutagenic GelRed DNA dye, for one-step DNA gel loading and staining. This is an improved version of our original 6X GelRed loading buffer that reduces the effects of dye on DNA migration, for bright si |
| Source                     | Synthetic  |