

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **CMV-p65(CMV100), CF405S conjugate, 0.1mg/mL, Clone: [CMV100], Mouse, Monoclonal BOT-BNC040100-100**

Article Name	CMV-p65(CMV100), CF405S conjugate, 0.1mg/mL, Clone: [CMV100], Mouse, Monoclonal
Biozol Catalog Number	BOT-BNC040100-100
Supplier Catalog Number	BNC040100-100
Alternative Catalog Number	BOT-BNC040100-100-100UL
Manufacturer	Biotium
Host	Mouse
Category	Antikörper
Species Reactivity	Virus
Immunogen	Recombinant CMV-p65 protein
Conjugation	CF405S
Product Description	Ten to forty percent of the patients with acquired immunodeficiency syndrome (AIDS) develop cytomegalovirus (CMV) infections. In some patients with AIDS, CMV is detected in the bronchoalveolar lavage fluid (BALF), urine, and other specimens, even whe...
Clonality	Monoclonal
Concentration	0.1 mg/mL
Clone Designation	[CMV100]
Molecular Weight	65 kDa

UniProt	Not Applicable
Buffer	PBS, 0.1% BSA, 0.05% azide
Source	Animal
Application Notes	For coating for ELISA, order Ab without BSA Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Optimal dilution and staining procedure for a specific application should be determined by user Recommended starting concentrations for titration are 1-2 ug/mL for most applications, or 1 ug/million cells/100 uL for flow cytometry