

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Melanoma Marker(HMB45 + M2-7C10 + M2-9E3 + T311), 0.2mg/mL, Clone: [HMB45 M2-7C10 M2-9E3 T311], Mouse, Monoclonal BOT-BNUB0701-100

Article Name	Melanoma Marker(HMB45 + M2-7C10 + M2-9E3 + T311), 0.2mg/mL, Clone: [HMB45 M2-7C10 M2-9E3 T311], Mouse, Monoclonal
Biozol Catalog Number	BOT-BNUB0701-100
Supplier Catalog Number	BNUB0701-100
Alternative Catalog Number	BOT-BNUB0701-100-100UL
Manufacturer	Biotium
Host	Mouse
Category	Antikörper
Application	IHC, WB
Species Reactivity	Human
Immunogen	Recombinant hMART-1 protein (M2-7C10, M2-9E3), Recombinant tyrosinase protein (T311), Extract of pigmented melanoma metastases from lymph nodes (HMB45)
Product Description	This antibody cocktail recognizes three melanoma-specific proteins, which include MART-1, Tyrosinase and gp100. MART-1 is a newly identified melanocyte differentiation antigen recognized by autologous cytotoxic T lymphocytes. Tyrosinase is one of the...
Clonality	Monoclonal
Concentration	0.2 mg/mL
Clone Designation	[HMB45 M2-7C10 M2-9E3 T311]

Molecular Weight	20-22 kDa (doublet) (MART), 70-80 kDa (Tyrosinase), 90-100 kDa (gp100)
UniProt	Q16655
Buffer	PBS, 0.05% BSA, 0.05% azide
Source	Animal
Application Notes	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistology (Formalin-fixed) 0.5-1.0 ug/mL Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user