

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

gp100 / Melanosome / PMEL17 / SILV (Melanoma Marker)(HMB45 + PMEL/783), 0.2mg/mL, Clone: [HMB45 PMEL/783], Mouse, Monoclonal BOT-BNUB0951-100

Article Name	gp100 / Melanosome / PMEL17 / SILV (Melanoma Marker)(HMB45 + PMEL/783), 0.2mg/mL, Clone: [HMB45 PMEL/783], Mouse, Monoclonal
Biozol Catalog Number	BOT-BNUB0951-100
Supplier Catalog Number	BNUB0951-100
Alternative Catalog Number	BOT-BNUB0951-100-100UL
Manufacturer	Biotium
Host	Mouse
Category	Antikörper
Application	IHC, WB
Species Reactivity	Human
Immunogen	Extract of pigmented melanoma metastases from lymph nodes (HMB45), Recombinant human SILV protein (PMEL/783)
Product Description	By immunohistochemistry, this antibody specifically recognizes a protein in melanocytes and melanomas. This MAb reacts with junctional and blue nevus cells and variably with fetal and neonatal melanocytes. Intradermal nevi, normal adult melanocytes, ...
Clonality	Monoclonal
Concentration	0.2 mg/mL
Clone Designation	[HMB45 PMEL/783]
Molecular Weight	90-100 kDa

UniProt	P40967
Buffer	PBS, 0.05% BSA, 0.05% azide
Source	Animal
Application Notes	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistochemistry (formalin-fixed): 0.5-1.0 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user