

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

**TIMP2 (Tissue Inhibitor of Metalloproteinase 2) (TIMP2/2044), 0.2mg/mL,  
Clone: [TIMP2/2044], Mouse, Monoclonal  
BOT-BNUB2044-100**

|                            |   |
|----------------------------|---|
| Article Name               | TIMP2 (Tissue Inhibitor of Metalloproteinase 2) (TIMP2/2044), 0.2mg/mL, Clone: [TIMP2/2044], Mouse, Monoclonal  |
| Biozol Catalog Number      | BOT-BNUB2044-100  |
| Supplier Catalog Number    | BNUB2044-100  |
| Alternative Catalog Number | BOT-BNUB2044-100-100UL  |
| Manufacturer               | Biotium   |
| Host                       | Mouse   |
| Category                   | Antikörper  |
| Application                | IHC   |
| Species Reactivity         | Human   |
| Immunogen                  | Recombinant human TIMP2 protein fragment (aa 48-220) (exact sequence is proprietary)  |
| Product Description        | This antibody recognizes a protein of 21 kDa, identified as tissue inhibitor of metalloproteinases-2 (TIMP-2). It is closely related to TIMP-1 and shows the highest binding affinity to both the latent (pro) and active forms of 72 kDa Type IV collagen... |
| Clonality                  | Monoclonal  |
| Concentration              | 0.2 mg/mL   |
| Clone Designation          | [TIMP2/2044]  |
| Molecular Weight           | 21 kDa  |

|                   |  |
|-------------------|--|
| UniProt           | P16035   |
| Buffer            | PBS, 0.05% BSA, 0.05% azide  |
| Source            | Animal   |
| Application Notes | Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody ELISA: 2-4 ug/m for coating order Ab without BSA Immunohistology (formalin): 1-2 ug/mL for 30 minutes at RT Staining of formalin-fixed tissues requires boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user |