

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Herpes Virus 8 (HHV8)(HHV8/3606), 1mg/mL, Clone: [HHV8/3606], Mouse, Monoclonal BOT-BNUM3606-50

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name | Human Herpes Virus 8 (HHV8)(HHV8/3606), 1mg/mL, Clone: [HHV8/3606], Mouse, Monoclonal |
| Biozol Catalog Number | BOT-BNUM3606-50 |
| Supplier Catalog Number | BNUM3606-50 |
| Alternative Catalog Number | BOT-BNUM3606-50-50UL |
| Manufacturer | Biotium |
| Host | Mouse |
| Category | Antikörper |
| Application | IHC |
| Species Reactivity | Virus |
| Immunogen | Recombinant protein corresponding to the latent nuclear antigen 1 molecule of HHV8 |
| Product Description | HHV 8 encodes a latent nuclear antigen (LNA), which is the product of the viral gene orf 73. LNA is capable of forming a complex with retinoblastoma susceptibility gene product, which may be related to its oncogenic activity. HHV8 is associated with ... |
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Clone Designation | [HHV8/3606] |
| Molecular Weight | Not Known |

| | |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| UniProt | Not Applicable |
| Buffer | PBS, no BSA, no azide |
| Source | Animal |
| Application Notes | ELISA: For coating use Ab at 1-5 ug/mL, order Ab without BSA, Optimal dilution for a specific application should be determined. Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody |