

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## DDDDK-tag polyclonal antibody, Unconjugated, Rabbit BWT-AP0007

Article Name	DDDDK-tag polyclonal antibody, Unconjugated, Rabbit
Biozol Catalog Number	BWT-AP0007
Supplier Catalog Number	AP0007
Alternative Catalog Number	BWT-AP0007-50UL,BWT-AP0007-100UL
Manufacturer	Bioworld Technology
Host	Rabbit
Category	Antikörper
Application	IF, IP, WB
Species Reactivity	AII
Immunogen	Synthetic peptide DDDDK-tag
Conjugation	Unconjugated
Product Description	Epitope tags are useful for the labeling and detection of proteins using immunoblotting, immunoprecipitation and immunostaining techniques. Due to their small size, they are unlikely to affect the tagged proteins biochemical properties
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Form	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Application Dilute	WB: 1:2000~1:5000 IP: 1:50~1:200 IF: 1:50~1:200

Application Notes

DDDDK-tag polyclonal antibody detects over-expressed or recombinant proteins containing the DDDDK epitope tag.