



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

Recombinant Human beta-Defensin 2 (rHuBD-2) BWT-PR1008

Recombinant Human beta-Defensin 2 (rHuBD-2)
BWT-PR1008
PR1008
BWT-PR1008-5UG,BWT-PR1008-20UG,BWT-PR1008-1.0MG
Bioworld Technology
Biochemikalien
Defensins (alpha and beta) are cationic peptides with a broad spectrum of antimicrobial activity that comprise an important arm of the innate immune system. The alpha-defensins are distinguished from the beta-defensins by the pairing of their three d
Approximately 4.3 KDa, a single non-glycosylated polypeptide chain containing 41 amino acids.
Escherichia coli.
>98% by SDS-PAGE and HPLC analyses.
Sterile Filtered White lyophilized (freeze-dried) powder.
GIGDPVTCLK SGAICHPVFC PRRYKQIGTC GLPGTKCCKK P
Lyophilized from a 0.2m filtered concentrated solution in 20mM PB, pH 7.4, 130mM NaCl.
We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportio