

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Sonic Hedgehog (Shh) Murine E. coli Preis auf Anfrage BVD-RBG20123100

Article Name	Sonic Hedgehog (Shh) Murine E. coli Preis auf Anfrage
Biozol Catalog Number	BVD-RBG20123100
Supplier Catalog Number	RBG20123100
Alternative Catalog Number	BVD-RBG20123100-100
Manufacturer	BioVendor
Category	Proteine/Peptide
Species Reactivity	Mouse
Product Description	Members of the Hedgehog (Hh) family are highly conserved proteins that are widely represented throughout the animal kingdom. The three known mammalian Hh proteins, Sonic (Shh), Desert (Dhh) and Indian (Ihh), are structurally related, and share a high
Source	E. coli
Expression System	E. coli
Purity	90%
Sequence	IVIGPGRGFG KRRHPKKLTP LAYKQFIPNV AEKTLGASGR YEGKITRNSE RFKELTPNYN PDIIFKDEEN TGADRLMTQR CKDKLNALAI SVMNQWPGVK LRVTEGWDED GHHSEESLHY EGRAVDITTS DRDRSKYGML ARLAVEAGFD WVYYESKAHI HCSVKAENSVAAKSGG
Storage	-20°C
Target	Sonic Hedgehog

## Application Dilute

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1–0.5 mg/ml. Do not vortex. For extended storage, it is recommended to further dilute in a buffer containing a carrier protein (example 0.1% BSA) and store in working...