

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## ELISA Kit for Myosin Light Chain 6, Alkali, Smooth Muscle And Non Muscle (MYL6) USC-SED426MU

Article Name	ELISA Kit for Myosin Light Chain 6, Alkali, Smooth Muscle And Non Muscle (MYL6)
Biozol Catalog Number	USC-SED426MU
Supplier Catalog Number	SED426Mu
Alternative Catalog Number	USC-SED426MU-48,USC-SED426MU-96
Manufacturer	Cloud-Clone
Category	Kits/Assays
Species Reactivity	Mouse
Range	0.156-10ng/mL
Sensitivity	The minimum detectable dose of this kit is typically less than 0.060ng/mL
Samples	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Application Notes	Enzyme-linked immunosorbent assay for Antigen Detection.