

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## RH IL-2c126s (AF) SHD-100-222AF-1MG

Article Name	RH IL-2c126s (AF)
Biozol Catalog Number	SHD-100-222AF-1MG
Supplier Catalog Number	100-222AF-1MG
Alternative Catalog Number	SHD-100-222AF-1MG
Manufacturer	Fujifilm Irvine Scientific
Category	Proteine/Peptide
Application	Bioassay, SDS-PAGE
Species Reactivity	Human
Product Description	IL-2C126S, an immune-regulating cytokine originating from lymphocytes, operates via the IL-2R receptor to trigger the expansion of activated T cells and encourage T cell maturation (1). IL-2 mediates its action by binding to IL-2 receptors (IL-2R), c
Molecular Weight	15.5
UniProt	P60568
Buffer	0.1% Trifluoroacetic Acid (TFA)
Source	E.coli
Purity	97
Form	Lyophilized from a 0.2 μm filtered solution

Sequence	MAPTSSSTKK TQLQLEHLLL DLQMILNGIN NYKNPKLTRM LTFKFYMPKK ATELKHLQCL EEELKPLEEV LNLAQSKNFH LRPRDLISNI NVIVLELKGS ETTFMCEYAD ETATIVEFLN RWITFSQSII STLT
Application Notes	Centrifuge vial before opening. Suspend the product by gently pipetting the above recommended solution down the sides of the vial. DO NOT VORTEX. Allow several minutes for complete reconstitution. For prolonged storage, dilute to working aliquots in a 0.