

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GENLISA Human Solute Carrier Family 2, Facilitated Glucose Transporter Member 1 (SLC2A1) ELISA KBS-KBH0319

Article Name	GENLISA Human Solute Carrier Family 2, Facilitated Glucose Transporter Member 1 (SLC2A1) ELISA
Biozol Catalog Number	KBS-KBH0319
Supplier Catalog Number	KBH0319
Alternative Catalog Number	KBS-KBH0319-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Human
Product Description	The GENLISA Human Solute Carrier Family 2, Facilitated Glucose Transporter Member 1 (SLC2A1) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Ac
Range	0.156 - 10 ng/ml
Sensitivity	0.05 ng/ml
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.
Application Notes	Colorimetric, Absorbance measured at 450nm