

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## GENLISA Human Zinc Transporter 3 / solute carrier family 30 member 3 (SLC30A3) ELISA KBS-KBH0370

Article Name	GENLISA Human Zinc Transporter 3 / solute carrier family 30 member 3 (SLC30A3) ELISA
Biozol Catalog Number	KBS-KBH0370
Supplier Catalog Number	KBH0370
Alternative Catalog Number	KBS-KBH0370-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Human
Product Description	The GENLISA Human Zinc Transporter 3 / solute carrier family 30 member 3 (SLC30A3) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and
Range	0.156 - 10 ng/ml
Sensitivity	0.05 ng/ml
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.
Application Notes	Colorimetric, Absorbance measured at 450nm