

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GENLISA Human G Protein-Activated Inward Rectifier Potassium Channel 2 (KCNJ6) ELISA KBS-KBH11182

Article Name	GENLISA Human G Protein-Activated Inward Rectifier Potassium Channel 2 (KCNJ6) ELISA
Biozol Catalog Number	KBS-KBH11182
Supplier Catalog Number	KBH11182
Alternative Catalog Number	KBS-KBH11182-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Human
Product Description	The GENLISA Human G Protein-Activated Inward Rectifier Potassium Channel 2 (KCNJ6) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and
Range	0 - 1000 ng/ml* (May vary. Please confirm)
Sensitivity	30 ng/ml* (May vary. Please confirm)
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.
Application Notes	Colorimetric, Absorbance measured at 450nm