

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Proteinase-Activated Receptor 4 (F2RL3) ELISA KBS-KBH1801

|                            |                                                                                                                                                                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Proteinase-Activated Receptor 4 (F2RL3) ELISA                                                                                                                                                                                                    |
| Biozol Catalog Number      | KBS-KBH1801                                                                                                                                                                                                                                                    |
| Supplier Catalog Number    | KBH1801                                                                                                                                                                                                                                                        |
| Alternative Catalog Number | KBS-KBH1801-1                                                                                                                                                                                                                                                  |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                            |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                    |
| Species Reactivity         | Human                                                                                                                                                                                                                                                          |
| Product Description        | The GENLISA Human Proteinase-Activated Receptor 4 (F2RL3) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precision Validated again... |
| Range                      | 0.781 - 50 ng/ml                                                                                                                                                                                                                                               |
| Sensitivity                | 0.35 ng/ml                                                                                                                                                                                                                                                     |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                             |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                     |