

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GENLISA Human Putative Sodium-Coupled Neutral Amino Acid Transporter 10 (SLC38A10) ELISA KBS-KBH6045

Article Name	GENLISA Human Putative Sodium-Coupled Neutral Amino Acid Transporter 10 (SLC38A10) ELISA
Biozol Catalog Number	KBS-KBH6045
Supplier Catalog Number	KBH6045
Alternative Catalog Number	KBS-KBH6045-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Human
Product Description	The GENLISA Putative Sodium-Coupled Neutral Amino Acid Transporter 10 (SLC38A10) ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistency - Accura...
Range	0 - 1000 ng/ml* (May vary. Please confirm)
Sensitivity	30 ng/ml* (May vary. Please confirm)
Samples	Serum, Plasma, Cell Culture Supernatant and other biological samples. Validated for citrated/EDTA Plasma only.
Application Notes	Colorimetric, Absorbance measured at 450nm