

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GENLISA Human Sodium-Dependent Serotonin Transporter (SLC6A4) ELISA KBS-KBH6189

Article Name	GENLISA Human Sodium-Dependent Serotonin Transporter (SLC6A4) ELISA
Biozol Catalog Number	KBS-KBH6189
Supplier Catalog Number	KBH6189
Alternative Catalog Number	KBS-KBH6189-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Human
Product Description	The GENLISA Sodium-Dependent Serotonin Transporter (SLC6A4) ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistency - Accuracy and Precision Vali...
Range	0 - 1000 ng/ml* (May vary. Please confirm)
Sensitivity	30 ng/ml* (May vary. Please confirm)
Samples	Serum, Plasma, Cell Culture Supernatant and other biological samples. Validated for citrated/EDTA Plasma only.
Application Notes	Colorimetric, Absorbance measured at 450nm