

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

GENLISA Human Anti-SARS-CoV-2 (Covid-19) IgG to spike S1+S2 Delta/Delta Plus Variant (B.1.617.2, AY.1,AY.2) Quantitative TITRATION ELISA KBS-KBVH015-46

| | |
|----------------------------|--|
| Article Name | GENLISA Human Anti-SARS-CoV-2 (Covid-19) IgG to spike S1+S2 Delta/Delta Plus Variant (B.1.617.2, AY.1,AY.2) Quantitative TITRATION ELISA |
| Biozol Catalog Number | KBS-KBVH015-46 |
| Supplier Catalog Number | KBVH015-46 |
| Alternative Catalog Number | KBS-KBVH015-46-1 |
| Manufacturer | Krishgen Biosystems |
| Category | Kits/Assays |
| Product Description | GENLISA Human Anti-SARS-CoV-2 (Covid-19) IgG to spike S1+S2 Delta/Delta Plus Variant (B.1.617.2, AY.1,AY.2) Quantitative TITRATION ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High... |
| Range | 0 - 1000 ng/ml |
| Sensitivity | 15 ng/ml |
| Samples | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further. Note the kit has been validated using human matrices only. For any other animal matrix, plea |
| Application Dilute | Dilute at 1:1000 with the Sample Diluent provided in the kit. Please refer to the IFU enclosed with the kit supplied for specific guidance. |