

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **GENLISA Mouse Epithelial Neutrophil Activating Peptide 78 (ENA78) ELISA KBS-KLM0427**

Article Name	GENLISA Mouse Epithelial Neutrophil Activating Peptide 78 (ENA78) ELISA
Biozol Catalog Number	KBS-KLM0427
Supplier Catalog Number	KLM0427
Alternative Catalog Number	KBS-KLM0427-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Mouse
Product Description	This GENLISA Mouse Epithelial Neutrophil Activating Peptide 78 (ENA78) ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistency - Accuracy and Pre...
Range	62.5 - 4,000 pg/ml
Sensitivity	27.2 pg/ml
Samples	Serum, plasma and tissue homogenates. For other matrices, please connect with our technical team