

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

PurKine(TM) MBP-Tag Dextrin Resin ASC-BMR2020

Article Name	PurKine(TM) MBP-Tag Dextrin Resin
Biozol Catalog Number	ASC-BMR2020
Supplier Catalog Number	BMR2020
Alternative Catalog Number	ASC-BMR2020-5, ASC-BMR2020-25
Manufacturer	Abbkine Scientific
Category	Sonstiges
Product Description	PurKine(TM) MBP-Tag Protein Purification Dextrin Resin 4FF
Tag	MBP Tag
Buffer	Liquid in PBS buffer containing 20% ethanol.
Form	Liquid solution with 50% slurry in 20% ethanol.
Application Notes	PurKine(TM) MBP-Tag Dextrin Resin effectively purifies high levels of overexpressed MBP-tagged fusion proteins in gravity column procedures at a variety of scales. The Resin consists of 90µm beads of highly cross-linked 4% agarose, to which Dextrin has been