

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

CheKine(TM) Micro Plant Anthocyanin Content Assay Kit ASC-KTB3011

Article Name	CheKine(TM) Micro Plant Anthocyanin Content Assay Kit
Biozol Catalog Number	ASC-KTB3011
Supplier Catalog Number	KTB3011
Alternative Catalog Number	ASC-KTB3011-48, ASC-KTB3011-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Product Description	CheKine(TM)CheKine(TM) Micro Plant Anthocyanin Content Assay Kit provides a simple, sensitive, rapid anthocyanin content detection method in plant tissues. Detailed sample preparation and result calculation methods are provided
Tag	Anthocyanin Content
Application Notes	It is recommended to use fresh samples, which can be stored at -80 C for several weeks if the experiment is not carried out immediately. When measuring, the temperature and time of thawing should be controlled. When thawing at room temperature, the sa