

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rat Tenascin-R (TNR) ELISA Kit ASC-KTE100053

Article Name	Rat Tenascin-R (TNR) ELISA Kit
Biozol Catalog Number	ASC-KTE100053
Supplier Catalog Number	KTE100053
Alternative Catalog Number	ASC-KTE100053-48,ASC-KTE100053-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Rat
Product Description	This Rat Tenascin-R (TNR) ELISA Kit employs a two-site sandwich ELISA to quantitate TNR....
Range	Please inquire
Sensitivity	Please inquire
Tag	TNR
NCBI	<a href="#">25567</a>
UniProt	<a href="#">Q05546</a>
Samples	Cell culture supernatants Serum Plasma Other biological fluids

Application Notes

This Rat Tenascin-R (TNR) ELISA Kit employs a two-site sandwich ELISA to quantitate TNR in samples. An antibody specific for TNR has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTNR present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNR is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNR bound in the initial step. The color development is stopped and the intensity of the color is measured.