

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Rat Trypsinogen activation peptide (TAP) ELISA Kit ASC-KTE100184

| Article Name               | Rat Trypsinogen activation peptide (TAP) ELISA Kit   |
|----------------------------|--|
| Biozol Catalog Number      | ASC-KTE100184  |
| Supplier Catalog Number    | KTE100184  |
| Alternative Catalog Number | ASC-KTE100184-48, ASC-KTE100184-96   |
| Manufacturer               | Abbkine Scientific   |
| Category                   | Kits/Assays  |
| Application                | ELISA  |
| Species Reactivity         | Rat  |
| Product Description        | This Rat Trypsinogen activation peptide (TAP) ELISA Kit employs a two-site sandwich ELISA to quantitate TAP  |
| Range                      | Please inquire   |
| Sensitivity                | Please inquire   |
| Tag                        | TAP  |
| Samples                    | Cell culture supernatantsSerumPlasmaOther biological fluids  |
| Application Notes          | This Rat Trypsinogen activation peptide (TAP) ELISA Kit employs a two-site sandwich ELISA to quantitate TAP in samples. An antibody specific for TAP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTAP prese |