

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse Transmembrane protein 17 (TMEM17) ELISA Kit ASC-KTE70183

Article Name	Mouse Transmembrane protein 17 (TMEM17) ELISA Kit
Biozol Catalog Number	ASC-KTE70183
Supplier Catalog Number	KTE70183
Alternative Catalog Number	ASC-KTE70183-48,ASC-KTE70183-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse Transmembrane protein 17 (TMEM17) ELISA Kit employs a two-site sandwich ELISA to quantitate TMEM17....
Range	Please inquire
Sensitivity	Please inquire
Tag	TMEM17
NCBI	<a href="#">103765</a>
UniProt	<a href="#">Q8K0U3</a>
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	<p>This Mouse Transmembrane protein 17 (TMEM17) ELISA Kit employs a two-site sandwich ELISA to quantitate TMEM17 in samples. An antibody specific for TMEM17 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TMEM17 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TMEM17 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TMEM17 bound in the initial step. The color development is stopped and the intensity of the color is measured.</p>
-------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------