

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Serine palmitoyltransferase 2 (SPTLC2) ELISA Kit ASC-KTE70251

Article Name	Mouse Serine palmitoyltransferase 2 (SPTLC2) ELISA Kit
Biozol Catalog Number	ASC-KTE70251
Supplier Catalog Number	KTE70251
Alternative Catalog Number	ASC-KTE70251-48,ASC-KTE70251-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse Serine palmitoyltransferase 2 (SPTLC2) ELISA Kit employs a two-site sandwich ELISA to quantitate SPTLC2
Range	Please inquire
Sensitivity	Please inquire
Tag	SPTLC2
NCBI	20773
UniProt	P97363
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes

This Mouse Serine palmitoyltransferase 2 (SPTLC2) ELISA Kit employs a two-site sandwich ELISA to quantitate SPTLC2 in samples. An antibody specific for SPTLC2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anySPTLC2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SPTLC2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SPTLC2 bound in the initial step. The color development is stopped and the intensity of the color is measured.