

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse Signal transducer and activator of transcription 3 (STAT3) ELISA Kit ASC-KTE70264

|                            |   |
|----------------------------|---|
| Article Name               | Mouse Signal transducer and activator of transcription 3 (STAT3) ELISA Kit  |
| Biozol Catalog Number      | ASC-KTE70264  |
| Supplier Catalog Number    | KTE70264  |
| Alternative Catalog Number | ASC-KTE70264-48,ASC-KTE70264-96   |
| Manufacturer               | Abbkine Scientific  |
| Category                   | Kits/Assays   |
| Application                | ELISA   |
| Species Reactivity         | Mouse   |
| Product Description        | This Mouse Signal transducer and activator of transcription 3 (STAT3) ELISA Kit employs a two-site sandwich ELISA to quantitate STAT3.... |
| Range                      | Please inquire  |
| Sensitivity                | Please inquire  |
| Tag                        | STAT3   |
| NCBI                       | 20848   |
| UniProt                    | P42227  |
| Samples                    | Cell culture supernatants<br>Serum<br>Plasma<br>Other biological fluids   |

Application Notes

This Mouse Signal transducer and activator of transcription 3 (STAT3) ELISA Kit employs a two-site sandwich ELISA to quantitate STAT3 in samples. An antibody specific for STAT3 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anySTAT3 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for STAT3 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of STAT3 bound in the initial step. The color development is stopped and the intensity of the color is measured.