

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## IPKine(TM) Anti-DDDDK Magnetic IP Kit, Magnetic Beads ASC-KTI2014

Article Name	IPKine(TM) Anti-DDDDK Magnetic IP Kit, Magnetic Beads
Biozol Catalog Number	ASC-KTI2014
Supplier Catalog Number	KTI2014
Alternative Catalog Number	ASC-KTI2014-20, ASC-KTI2014-100
Manufacturer	Abbkine Scientific
Category	Antikörper
Application	IP
Species Reactivity	Bacteria, Mammal
Immunogen	Synthetic Peptide
Conjugation	Magnetic Beads
Product Description	The kit contains high-quality Anti-Flag magnetic beads and the necessary reagents for immunoprecipitation that have been optimized and validated, making immunoprecipitation (IP) experiments simpler, more convenient, and more efficient
Tag	Flag
Application Notes	Anti-DDDDK-Tag (Flag-Tag) Magnetic Beads are prepared by covalently coupling Anti-Flag-Tag Mouse Monoclonal Antibody to crosslinked Magnetic Beads. The Optimized Anti-Flag Magnetic Beads have more efficient antigen binding capacity. This kit provides thr