



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

TDP-43 WT BacMam - 30mL MTM-N0200N-030

Article Name	TDP-43 WT BacMam - 30mL
Biozol Catalog Number	MTM-N0200N-030
Supplier Catalog Number	N0200N-030
Alternative Catalog Number	MTM-N0200N-030
Manufacturer	Montana Molecular
Category	Molekularbiologie
Product Description	Wild-type TDP-43 packaged in BacMam for expression in mammalian cells. BacMam provides even, near 100% transduction efficiency in HEK293 cells as well as expression in many other cell lines, iPSCs, and primary cultures. Expression is easily titrateab