

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Brilliant blue G-250, CAS [[6104-58-1]] MCE-HY-D0014

Article Name	Brilliant blue G-250, CAS [[6104-58-1]]
Biozol Catalog Number	MCE-HY-D0014
Supplier Catalog Number	HY-D0014
Alternative Catalog Number	MCE-HY-D0014-100MG,MCE-HY-D0014-500MG,MCE-HY-D0014- 250MG,MCE-HY-D0014-1G
Manufacturer	MedchemExpress
Category	Kits/Assays
Product Description	Brilliant Blue G-250 is a dye commonly used for the visualization of proteins separated by SDS-PAGE, offering a simple staining procedure and high quantitation. In the Bradford protein assay, protein concentrations are determined by the absorbance at
Molecular Weight	854.02
Purity	95
CAS Number	[6104-58-1]
Formula	C47H48N3NaO7S2
Target	Biochemical Assay Reagents
Application Notes	MCE Product type: Biochemical Assay Reagents