

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

DPY19L2 Human Pre-designed siRNA Set A MCE-HY-RS03976

Article Name	DPY19L2 Human Pre-designed siRNA Set A
Biozol Catalog Number	MCE-HY-RS03976
Supplier Catalog Number	HY-RS03976
Alternative Catalog Number	MCE-HY-RS03976-1SET
Manufacturer	MedchemExpress
Category	Molekularbiologie
Product Description	DPY19L2 Human Pre-designed siRNA Set A contains three designed siRNAs for DPY19L2 gene (Human), as well as a negative control, a positive control, and a FAM-labeled negative control. Components DPY19L2 siRNA-1: 5 nmol (HPLC) DPY19L2 siRNA-2: 5 nmol (HPLC) DPY19L2 siRNA-3: 5 nmol (HPLC) DPY19L2 siRNA-4: 5 nmol (HPLC) DPY19L2 siRNA-5: 5 nmol (HPLC) DPY19L2 siRNA-6: 5 nmol (HPLC) DPY19L2 siRNA-7: 5 nmol (HPLC) DPY19L2 siRNA-8: 5 nmol (HPLC) DPY19L2 siRNA-9: 5 nmol (HPLC) DPY19L2 siRNA-10: 5 nmol (HPLC) DPY19L2 siRNA-11: 5 nmol (HPLC) DPY19L2 siRNA-12: 5 nmol (HPLC) DPY19L2 siRNA-13: 5 nmol (HPLC) DPY19L2 siRNA-14: 5 nmol (HPLC) DPY19L2 siRNA-15: 5 nmol (HPLC) DPY19L2 siRNA-16: 5 nmol (HPLC) DPY19L2 siRNA-17: 5 nmol (HPLC) DPY19L2 siRNA-18: 5 nmol (HPLC) DPY19L2 siRNA-19: 5 nmol (HPLC) DPY19L2 siRNA-20: 5 nmol (HPLC) DPY19L2 siRNA-21: 5 nmol (HPLC) DPY19L2 siRNA-22: 5 nmol (HPLC) DPY19L2 siRNA-23: 5 nmol (HPLC) DPY19L2 siRNA-24: 5 nmol (HPLC) DPY19L2 siRNA-25: 5 nmol (HPLC) DPY19L2 siRNA-26: 5 nmol (HPLC) DPY19L2 siRNA-27: 5 nmol (HPLC) DPY19L2 siRNA-28: 5 nmol (HPLC) DPY19L2 siRNA-29: 5 nmol (HPLC) DPY19L2 siRNA-30: 5 nmol (HPLC) DPY19L2 siRNA-31: 5 nmol (HPLC) DPY19L2 siRNA-32: 5 nmol (HPLC) DPY19L2 siRNA-33: 5 nmol (HPLC) DPY19L2 siRNA-34: 5 nmol (HPLC) DPY19L2 siRNA-35: 5 nmol (HPLC) DPY19L2 siRNA-36: 5 nmol (HPLC) DPY19L2 siRNA-37: 5 nmol (HPLC) DPY19L2 siRNA-38: 5 nmol (HPLC) DPY19L2 siRNA-39: 5 nmol (HPLC) DPY19L2 siRNA-40: 5 nmol (HPLC) DPY19L2 siRNA-41: 5 nmol (HPLC) DPY19L2 siRNA-42: 5 nmol (HPLC) DPY19L2 siRNA-43: 5 nmol (HPLC) DPY19L2 siRNA-44: 5 nmol (HPLC) DPY19L2 siRNA-45: 5 nmol (HPLC) DPY19L2 siRNA-46: 5 nmol (HPLC) DPY19L2 siRNA-47: 5 nmol (HPLC) DPY19L2 siRNA-48: 5 nmol (HPLC) DPY19L2 siRNA-49: 5 nmol (HPLC) DPY19L2 siRNA-50: 5 nmol (HPLC) DPY19L2 siRNA-51: 5 nmol (HPLC) DPY19L2 siRNA-52: 5 nmol (HPLC) DPY19L2 siRNA-53: 5 nmol (HPLC) DPY19L2 siRNA-54: 5 nmol (HPLC) DPY19L2 siRNA-55: 5 nmol (HPLC) DPY19L2 siRNA-56: 5 nmol (HPLC) DPY19L2 siRNA-57: 5 nmol (HPLC) DPY19L2 siRNA-58: 5 nmol (HPLC) DPY19L2 siRNA-59: 5 nmol (HPLC) DPY19L2 siRNA-60: 5 nmol (HPLC) DPY19L2 siRNA-61: 5 nmol (HPLC) DPY19L2 siRNA-62: 5 nmol (HPLC) DPY19L2 siRNA-63: 5 nmol (HPLC) DPY19L2 siRNA-64: 5 nmol (HPLC) DPY19L2 siRNA-65: 5 nmol (HPLC) DPY19L2 siRNA-66: 5 nmol (HPLC) DPY19L2 siRNA-67: 5 nmol (HPLC) DPY19L2 siRNA-68: 5 nmol (HPLC) DPY19L2 siRNA-69: 5 nmol (HPLC) DPY19L2 siRNA-70: 5 nmol (HPLC) DPY19L2 siRNA-71: 5 nmol (HPLC) DPY19L2 siRNA-72: 5 nmol (HPLC) DPY19L2 siRNA-73: 5 nmol (HPLC) DPY19L2 siRNA-74: 5 nmol (HPLC) DPY19L2 siRNA-75: 5 nmol (HPLC) DPY19L2 siRNA-76: 5 nmol (HPLC) DPY19L2 siRNA-77: 5 nmol (HPLC) DPY19L2 siRNA-78: 5 nmol (HPLC) DPY19L2 siRNA-79: 5 nmol (HPLC) DPY19L2 siRNA-80: 5 nmol (HPLC) DPY19L2 siRNA-81: 5 nmol (HPLC) DPY19L2 siRNA-82: 5 nmol (HPLC) DPY19L2 siRNA-83: 5 nmol (HPLC) DPY19L2 siRNA-84: 5 nmol (HPLC) DPY19L2 siRNA-85: 5 nmol (HPLC) DPY19L2 siRNA-86: 5 nmol (HPLC) DPY19L2 siRNA-87: 5 nmol (HPLC) DPY19L2 siRNA-88: 5 nmol (HPLC) DPY19L2 siRNA-89: 5 nmol (HPLC) DPY19L2 siRNA-90: 5 nmol (HPLC) DPY19L2 siRNA-91: 5 nmol (HPLC) DPY19L2 siRNA-92: 5 nmol (HPLC) DPY19L2 siRNA-93: 5 nmol (HPLC) DPY19L2 siRNA-94: 5 nmol (HPLC) DPY19L2 siRNA-95: 5 nmol (HPLC) DPY19L2 siRNA-96: 5 nmol (HPLC) DPY19L2 siRNA-97: 5 nmol (HPLC) DPY19L2 siRNA-98: 5 nmol (HPLC) DPY19L2 siRNA-99: 5 nmol (HPLC) DPY19L2 siRNA-100: 5 nmol (HPLC)
Target	Small Interfering RNA (siRNA)
Application Notes	MCE Product type: Oligonucleotides