

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human TAP (Trypsinogen Activation Peptide) ELISA Kit EKL-EH24409

|                            |   |
|----------------------------|---|
| Article Name               | Human TAP (Trypsinogen Activation Peptide) ELISA Kit  |
| Biozol Catalog Number      | EKL-EH24409   |
| Supplier Catalog Number    | EH24409   |
| Alternative Catalog Number | EKL-EH24409-96T   |
| Manufacturer               | EnkiLife  |
| Category                   | Kits/Assays   |
| Species Reactivity         | Human   |
| Product Description        | This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Human TAP. During the reaction, Human TAP in the sample or standard competes with a fixed amount of Human TAP on the solid phase ... |
| Range                      | 0.124-10ng/mL   |
| Sensitivity                | 0.054ng/mL  |
| Samples                    | Serum,plasma and other biological fluids  |
| Target                     | TAP   |