

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Anti-DYKDDDDK-tag FAST IP kit ELA-EA-IP-K001

|                            |   |
|----------------------------|---|
| Article Name               | Anti-DYKDDDDK-tag FAST IP kit   |
| Biozol Catalog Number      | ELA-EA-IP-K001  |
| Supplier Catalog Number    | EA-IP-K001  |
| Alternative Catalog Number | ELA-EA-IP-K001-50T  |
| Manufacturer               | Elabscience   |
| Category                   | Kits/Assays   |
| Product Description        | Anti-DYKDDDDK-tag fast IP kit is composed of Anti-DYKDDDDK affinity gel, Mouse IgG affinity gel, DYKDDDDK-Tag Rabbit pAb and Goat Anti-Rabbit IgG(peroxidase/HRP conjugated), which is used for fast, efficient and specific (Co)IP of DYKDDDDK tag fusion... |