

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Cleaved PARP Mouse Monoclonal Antibody, Clone: [173] MBS8501787-0.1-CY7**

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Cleaved PARP Mouse Monoclonal Antibody, Clone: [173]                                                                                    |
| Biozol Catalog Number      | MBS8501787-0.1-CY7                                                                                                                      |
| Supplier Catalog Number    | MBS8501787-0.1-Cy7                                                                                                                      |
| Alternative Catalog Number | MBS8501787-0.1-CY7                                                                                                                      |
| Manufacturer               | MyBiosource                                                                                                                             |
| Category                   | Antikörper                                                                                                                              |
| Clonality                  | Monoclonal                                                                                                                              |
| Clone Designation          | [173]                                                                                                                                   |
| Form                       | In phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |