

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ToxinEraser(TM) Endotoxin Removal Kit GSC-L00338

Article Name	ToxinEraser(TM) Endotoxin Removal Kit
Biozol Catalog Number	GSC-L00338
Supplier Catalog Number	L00338
Alternative Catalog Number	GSC-L00338
Manufacturer	Genscript
Category	Proteine/Peptide
Storage	For storage of the column, wash the column with 10 ml of equilibration buffer and allow the column to drain completely. Add 1.5 ml of regeneration buffer supplemented with 0.02% sodium azide. Store at 2°C to 8°C. Do not freeze.