

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## AKIP Polyclonal Antibody, PerCP Conjugated, Rabbit BSS-BS-4287R-PERCP

AKIP Polyclonal Antibody, PerCP Conjugated, Rabbit
BSS-BS-4287R-PERCP
bs-4287R-PerCP
BSS-BS-4287R-PERCP-100
Bioss
Rabbit
Antikörper
IF
PerCP
May act as a negative regulator of Aurora-A kinase, by down-regulation through proteasome-dependent degradation
Polyclonal
1ug/ul
54998
Q9NWT8
Aqueous buffered solution containing 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
KLH conjugated synthetic peptide derived from human AKIP
Purified by Protein A.

Target	AKIP
Application Dilute	IF(IHC-P)(1:50-200), IF(IHC-F)(1:50-200), IF(ICC)(1:50-200)