

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Canine GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit NBS-EKX-B2QA1E-96

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| Article Name               | Canine GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor)<br>ELISA Kit    |
| Biozol Catalog Number      | NBS-EKX-B2QA1E-96                                                                |
| Supplier Catalog Number    | EKX-B2QA1E-96                                                                    |
| Alternative Catalog Number | NBS-EKX-B2QA1E-96                                                                |
| Manufacturer               | Nordic BioSite                                                                   |
| Category                   | Kits/Assays                                                                      |
| Application                | ELISA                                                                            |
| Species Reactivity         | Canine                                                                           |
| Product Description        | Canine GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor)<br>ELISA Kit... |
| Range                      | 15.625-1000pg/ml                                                                 |
| Sensitivity                | 9.375pg/ml                                                                       |
| NCBI                       | <a href="#">403923</a>                                                           |
| UniProt                    | <a href="#">P48749</a>                                                           |
| Target                     | GM-CSF/CSF2                                                                      |