

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Native Pseudomonas sp. Lipoprotein lipase, >20 U/mg, CAS [[08/02/9004]] CRM-DIA-210

Article Name	Native Pseudomonas sp. Lipoprotein lipase, >20 U/mg, CAS [[08/02/9004]]
Biozol Catalog Number	CRM-DIA-210
Supplier Catalog Number	DIA-210
Alternative Catalog Number	CRM-DIA-210-1
Manufacturer	Creative BioMart
Category	Biochemikalien
Concentration	GradeIII 20U/mg-solid or more (containing approx. 80% of stabilizers)
CAS Number	[08/02/9004]
Application Notes	This enzyme is useful for enzymatic determination of triglyceride in serum when coupled with L-alpha-glycerophosphate oxidase and glycerol kinase. Usually, the reaction can be completed in 5 minutes at 37C by using $2.5\sim3.0$ units of the enzyme per test (3.0ml) at pH around 7.0.