

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Insect Cathepsin(Cath) ELISA Kit CRM-IEK003

| | |
|----------------------------|---|
| Article Name | Insect Cathepsin(Cath) ELISA Kit |
| Biozol Catalog Number | CRM-IEK003 |
| Supplier Catalog Number | IEK003 |
| Alternative Catalog Number | CRM-IEK003-48, CRM-IEK003-96 |
| Manufacturer | Creative BioMart |
| Category | Kits/Assays |
| Application Notes | This kit adopts double antibody sandwich ELISA method. Biotinylated antibodies and horseradish peroxidase-labeled affinities are added sequentially. The TMB appears blue and turns yellow after the addition of termination solution. The OD value was measured at 450 nm, and the concentration of antigen was proportional to the OD450 value. |