



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Insect Membrane Protein, Palmitoylated 6(MPP6) ELISA Kit CRM-IEK008

Article Name	Insect Membrane Protein, Palmitoylated 6(MPP6) ELISA Kit
Biozol Catalog Number	CRM-IEK008
Supplier Catalog Number	IEK008
Alternative Catalog Number	CRM-IEK008-48, CRM-IEK008-96
Manufacturer	Creative BioMart
Category	Kits/Assays
Application Notes	This kit adopts double antibody sandwich ELISA method. Biotinylated antibodies and horseradish peroxidase-labeled affinities are added sequentially. The TMB appears blue and turns yellow after the addition of termination solution. The OD value was measured at 450 nm, and the concentration of antigen was proportional to the OD450 value.