

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Insect Acetoacetic Acid(ACAC) ELISA Kit CRM-IEK030

Article Name	Insect Acetoacetic Acid(ACAC) ELISA Kit
Biozol Catalog Number	CRM-IEK030
Supplier Catalog Number	IEK030
Alternative Catalog Number	CRM-IEK030-48, CRM-IEK030-96
Manufacturer	Creative BioMart
Category	Kits/Assays
Application Notes	This kit adopts double antibody sandwich ELISA method. Biotinylated antibodies and horseradish peroxidase-labeled affinities are added sequentially. The TMB appears blue and turns yellow after the addition of termination solution. The OD value was measured at 450 nm, and the concentration of antigen was proportional to the OD450 value.