

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Bacteroides fragilis N-acetylmuramic acid 6-phosphate etherase(murQ),Escherichia Coli Preis auf Anfrage CSB-EP683497BAAB

Article Name	Recombinant Bacteroides fragilis N-acetylmuramic acid 6-phosphate etherase(murQ),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP683497BAAB
Supplier Catalog Number	CSB-EP683497BAAB
Alternative Catalog Number	CSB-EP683497BAAB-1, CSB-EP683497BAAB-100, CSB-EP683497BAAB-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q5LI89
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MNSNIEKSDK PSFIKISEQP SLYDDLEKKS VREILEDINK EDQKVAIAVQ KAIPQIEKLV TQIVPRMKQG GRIFYMGAGT SGRLGVLDAS EIPPTFGMPP TLIIGLIAGG DTALRNPVEN AEDNTTRGWE ELTEHNINDK DTVIGIAASG TTPYVIGAMH AAREHGILTG CITSNPNSPM AAEDIPIEM IVGPEYVTGS SRMKSGTGQK MILNMITTSV MIQLGRVKGN KM

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.