

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Tachyglossus aculeatus aculeatus Amelogenin(AMEL),Escherichia Coli Preis auf Anfrage CSB-EP001664TKU

Article Name	Recombinant Tachyglossus aculeatus aculeatus Amelogenin(AMEL),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP001664TKU
Supplier Catalog Number	CSB-EP001664TKU
Alternative Catalog Number	CSB-EP001664TKU-1,CSB-EP001664TKU-100,CSB-EP001664TKU-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	O97647
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	LHHQIIPVLS QHQTPTHALQ SHHHIPVMAT QQPTQPPQPM MPMPGQHSVT PTQHHQSNLP QPGQQPFQPQ FPQKPTHRPI QPQAPVHPMP PMPQPQLPPM FPLQPLPPLL PDLPLEPWPA S

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.