

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Nitratiruptor sp. Ribosomal RNA large subunit methyltransferase N(rlmN), Escherichia Coli Preis auf Anfrage CSB-EP412354NFA

|                            |                                                                                                                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Nitratiruptor sp. Ribosomal RNA large subunit methyltransferase N(rlmN), Escherichia Coli Preis auf Anfrage                                                                                                                                                        |
| Biozol Catalog Number      | CSB-EP412354NFA                                                                                                                                                                                                                                                                |
| Supplier Catalog Number    | CSB-EP412354NFA                                                                                                                                                                                                                                                                |
| Alternative Catalog Number | CSB-EP412354NFA-1, CSB-EP412354NFA-100, CSB-EP412354NFA-20                                                                                                                                                                                                                     |
| Manufacturer               | Cusabio                                                                                                                                                                                                                                                                        |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                               |
| Species Reactivity         | Bacteria                                                                                                                                                                                                                                                                       |
| Tag                        | The tag type will be determined during production process.                                                                                                                                                                                                                     |
| UniProt                    | <a href="#">A6Q115</a>                                                                                                                                                                                                                                                         |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                    |
| Source                     | E.coli                                                                                                                                                                                                                                                                         |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                    |
| Form                       | Lyophilized powder                                                                                                                                                                                                                                                             |
| Sequence                   | MKNILDLTKE ELQQEVTPKF RANQIYQWIY QKGAKDFESM SNLPKSMREE<br>LKEKFTITPP KILNVEVSKD GSKKYLLGLQ DGHTVESVLL PMKKEERDEK<br>GNILKEARYT VCVSSQVGCK VGCEFCCLTAK GGFVRNLTPG EIVEQVLTIR<br>EDNNIPANRR VNIVYMGMGGE PLDNLENVAK AVKIFSDEHG MSISPRRQTI<br>STSG LAPKIK KLGEMNLGVL LAISLHAVDD EL |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|