

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Nitratiruptor sp. 30S ribosomal protein S17(rpsQ),Escherichia Coli Preis auf Anfrage CSB-EP412357NFA

Article Name	Recombinant Nitratiruptor sp. 30S ribosomal protein S17(rpsQ),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP412357NFA
Supplier Catalog Number	CSB-EP412357NFA
Alternative Catalog Number	CSB-EP412357NFA-1,CSB-EP412357NFA-100,CSB-EP412357NFA-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6Q1I7
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MAYKRVIQGK VVKKSGDKTV SLLVERKVVH PKYHKIVKRF KKYLVHDEHN KAKVGDIVTA VECRPISKRK SFRLKEIVQA GVE

Application Notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.