

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Recombinant Nitratiruptor sp. 50S ribosomal protein L28(rpmB),Escherichia Coli Preis auf Anfrage CSB-EP412360NFA

|                            |                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------|
| Article Name               | Recombinant Nitratiruptor sp. 50S ribosomal protein<br>L28(rpmB),Escherichia Coli Preis auf Anfrage |
| Biozol Catalog Number      | CSB-EP412360NFA                                                                                     |
| Supplier Catalog Number    | CSB-EP412360NFA                                                                                     |
| Alternative Catalog Number | CSB-EP412360NFA-1,CSB-EP412360NFA-100,CSB-EP412360NFA-20                                            |
| Manufacturer               | Cusabio                                                                                             |
| Category                   | Proteine/Peptide                                                                                    |
| Species Reactivity         | Bacteria                                                                                            |
| Tag                        | The tag type will be determined during production process.                                          |
| UniProt                    | <a href="#">A6Q1S4</a>                                                                              |
| Buffer                     | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                         |
| Source                     | E.coli                                                                                              |
| Purity                     | Greater than 85% as determined by SDS-PAGE.                                                         |
| Form                       | Lyophilized powder                                                                                  |
| Sequence                   | MSKKCQVTGK KPITGHNVSH ANNKTKRRFL PNIRTVRVTL PDGTKKKIKV<br>SAKGLRILKK KNFQVDL                        |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|