

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Recombinant Nitratiruptor sp. Exodeoxyribonuclease 7 large subunit(xseA),Escherichia Coli Preis auf Anfrage CSB-EP412363NFA

Article Name	Recombinant Nitratiruptor sp. Exodeoxyribonuclease 7 large subunit(xseA),Escherichia Coli Preis auf Anfrage
Biozol Catalog Number	CSB-EP412363NFA
Supplier Catalog Number	CSB-EP412363NFA
Alternative Catalog Number	CSB-EP412363NFA-1,CSB-EP412363NFA-100,CSB-EP412363NFA-20
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6Q232
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Source	E.coli
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Lyophilized powder
Sequence	MNVYTVSQLN EQIKNLLESH FVEVYVEGEV SRPTYHTSGH LYFSLKDEKS VIRCVMF RSA LAKVPFRVED GQKLIVAGKI GVIKPRGEYQ LYATELHP SG VGS LQLAFEQ LKAKLEKKG Y FASELKKPLP DFIQTIALVT SQTGAALQDM LR IIQNRWPL VKVYVVDTLV QGSDAAPMIA RSIAYADGLG VDVIVVGRGG GSLEDLWPFN EEIVADAIFE AKTPIVSAVG HE

Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-------------------	--